Transcend LCD Programming Base Station Patient Manual

ENGLISH

FRANÇAIS

DEUTSCH

ITALIANO

ESPAÑOL

NEDERLANDS

PORTUGUÊS

Description of Transcend LCD Programming Base Station

The Transcend LCD Programming Base Station (LCD) is a unique accessory to be used along with the Transcend Auto™, Transcend EZEX™, and Transcend® Sleep Apnea Therapy Starter System (PAP device). The LCD seats the PAP securely in the stand. A keypad with three buttons facilitates control and the screen visually displays compliance data. The LCD is powered by Somnetics-approved power supplies. The LCD and Patient manual are included in the package contents. Report any missing items to your homecare provider.

The LCD displays the following information to the patient:

Transcend Auto	Transcend EZEX	Transcend
<ul style="list-style-type: none"> Minimum Pressure in cm H2O 	<ul style="list-style-type: none"> Therapy Pressure in cm H2O 	<ul style="list-style-type: none"> Therapy Pressure in cm H2O
<ul style="list-style-type: none"> Maximum in cm H2O 	<ul style="list-style-type: none"> Ramp Duration in minutes 	<ul style="list-style-type: none"> Ramp Duration in minutes
<ul style="list-style-type: none"> Starting Pressure in cm H2O 	<ul style="list-style-type: none"> Ramp Pressure in cm H2O 	<ul style="list-style-type: none"> Ramp Pressure in cm H2O
<ul style="list-style-type: none"> Ramp Duration in minutes 	<ul style="list-style-type: none"> EZEX Pressure Relief 	<ul style="list-style-type: none"> Apnea Hypopnea Index (AHI)
<ul style="list-style-type: none"> Ramp Pressure in cm H2O 	<ul style="list-style-type: none"> Apnea Hypopnea Index (AHI) 	<ul style="list-style-type: none"> Leak in L/min
<ul style="list-style-type: none"> EZEX Pressure Relief 	<ul style="list-style-type: none"> Leak in L/min 	<ul style="list-style-type: none"> Device hours in both "Total Blower Hours" and "Total Patient Hours"
<ul style="list-style-type: none"> Apnea Hypopnea Index (AHI) 	<ul style="list-style-type: none"> Device hours in both "Total Blower Hours" and "Total Patient Hours" 	<ul style="list-style-type: none"> Device fault screen
<ul style="list-style-type: none"> Leak in L/min 	<ul style="list-style-type: none"> Device fault screen 	<ul style="list-style-type: none"> Save Screen
<ul style="list-style-type: none"> Device hours in both "Total Blower Hours" and "Total Patient Hours" 	<ul style="list-style-type: none"> Save Screen 	
<ul style="list-style-type: none"> Device fault screen 		
<ul style="list-style-type: none"> Save Screen 		

Note: The LCD will display a 15 second pressure average when in therapy mode. See "Therapy Mode" section for more information.

Note: Total Blower Hours and Total Patient Hours are divided into two separate screens.

Precautions for use

Warnings

- Use only those accessories provided or recommended by Somnetics. Use of accessories not recommended may result in increased electromagnetic emissions or decreased electromagnetic immunity of the unit.
- This device is not waterproof and therefore should not be exposed to environmental conditions where it may get wet.

Cautions

- Use of the LCD with Transcend P8 or P4 batteries may impact PAP therapy duration if only using battery power

Symbols

Consult instructions for use

Upper and lower temperature limits

Separate collection of electronic equipment per EC Directive 2002/96/EC - Waste Electrical and Electronic Equipment (WEEE)

Accessory should be kept dry

Attention: consult accompanying documents

Upper and lower humidity limits

Symbols and Abbreviations on LCD Display

Ramp

Hours

T

Total

Total Blower Hours on the PAP

Total Patient Hours on PAP

cmH₂O

PAP Pressure in cm H₂O

ℓ/min



Leak in Liters per minute

AHI

Apnea Hypopnea Index (AHI)

Setting up the LCD Programming Base Station

- Begin by connecting a Somnetics-approved power source into the LCD (Fig 1.) If using the AC power supply, plug it into the wall outlet.
- Place the PAP into the LCD (Fig 2.) by first sliding the front hooks of the PAP into the front slats of the LCD then press down on the back of the PAP until it snaps into place on the LCD.
- Connect the LCD to the PAP (Fig 3.) by using the USB cable and PAP power cable located at the back of the LCD.

Section 1: Operation

Patient Mode

1. The LCD must be connected to a power supply (e.g. AC power supply, Mobile Power Adaptor, or P4 or P8 Battery) for correct information to display and to change PAP ramp settings (see below). Follow the directions in 'Setting up the LCD Programming Base Station' section prior to operation.
2. To review therapy data, push the "→" button to cycle through the following screens in the order presented below. The "↑" or "↓" buttons do not change the display except as noted.
 - To turn on the LCD press any button.
 - To turn off the LCD, disconnect the power supply from the LCD. Remove the USB cable and PAP power cable from the PAP device.

Note: This device does not have a back button. If a screen is passed, the user must advance through the screens using the "→" button to return to the desired screen.

Note: Use the device in adequate lighting conditions in order to view the screens below.

Auto	EZEX	Transcend	Screen Display	Description
X			1  cmH ₂ O	Minimum Pressure (4-20 cm H ₂ O) 1 = min pressure; digits to the right display value
X			2  cmH ₂ O	Starting Pressure (min pressure - 20) 2 = max pressure; digits to the right display value
X			3  cmH ₂ O	Starting Pressure (between min and max pressure) 3 = starting pressure; digits to the right display value
	X	X	 cmH ₂ O	Therapy Pressure (4 – 20 cm H ₂ O)

Auto	EZEX	Transcend	Screen Display	Description
X	X	X		Ramp Duration (0-45 minutes) *Patient adjustable with “↑” or “↓” buttons (5 minute increments)
X	X	X	 cmH ₂ O	Ramp Pressure (4 cm H₂O – therapy/starting pressure) *Patient adjustable with “↑” or “↓” buttons (1 cm H ₂ O increments)
X	X		E2 	EZEX Pressure Relief Setting of 0, 1, 2, 3
X	X	X	 AHI	Apnea Hypopnea Index (AHI) From previous night *Initially displayed with 5 dashes until “↑” or “↓” buttons are pressed. Press either “↑” or “↓” to retrieve AHI information. Display will show “01234” until AHI value is retrieved.
X	X	X	 l/min	Leak from previous night *Initially displayed with 5 dashes until “↑” or “↓” buttons are pressed. Press either “↑” or “↓” to retrieve leak detection information. Display will show “01234” until AHI value is retrieved.
X	X	X	 hr →	Total Patient Hours /p
X	X	X	 hr → T	Total Blower Hours on Device
X	X	X	F 	Device Fault Code *Refer to Transcend PAP User Manuals for PAP Fault Codes
X	X	X	SAVE	Save

Therapy Mode

- Ensure that the hose is connected to the PAP and that your mask is attached. To prevent the PAP from shutting down you must cover the mask opening or put on the mask while the device is in operation.

Note: If there are any leaks the system will shut down.

- Press the power button on the PAP to start therapy mode.

	<p>Therapy Mode * When the PAP is in therapy mode (i.e. positive pressure delivery) this screen will continually provide a 15-second time average of the pressure sensed by the PAP unit</p>
---	---

Section 2: Power and Fault Conditions

1. The LCD must be fully connected as described in the 'Setting up the LCD Programming Base Station' section. If set up as described below, the LCD will continue to display previously programmed settings, but the information may not be correct and you will not be able to modify ramp settings.
 - Power supply is connected to the LCD and the PAP power cable is not connected to the PAP.
2. If the USB cable is not connected to a PAP device, the LCD will not display any settings.
3. When in Therapy mode, the screen will go blank after 45 seconds of inactivity.
4. Press any of the three buttons to return the display to the therapy pressure screen (Patient Mode) or the time average therapy pressure screen (Therapy Mode).
5. When the PAP is in a fault condition, the Therapy Mode is inaccessible. Manual reset of the PAP is required to obtain real-time pressure readings. (Refer to Transcend PAP user manuals for more information on fault conditions and manual reset of the device.)

Caring for your Transcend LCD Programming Base Station

Clean the exterior of the base station as needed. Use a damp cloth or paper towel with mild dishwashing soap to wipe clean. Never submerge in water. Do not get soap or water in base station.

Product Specifications

Dimensions (LxWxH)	8.0 in x 5.0 in x 3.25 in (20.3 cm x 12.7 cm x 8.3 cm)
Weight	0.72 lbs (326.0 g)
Operating Range	41°-95°F (5°-35°C) 10%-80% relative humidity (non-condensing) 0-8,000 ft (0-2,438 m)
Storage Range	-4° to 140°F (-20° to 60°C) 10%-90% relative humidity (non-condensing)

Manufacturer's declaration

Electromagnetic emissions

The LCD is intended for use in the electromagnetic environment specified below. The customer or the user of the system should ensure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment—guidance
RF radiated emissions CISPR 11	Group 1	The LCD uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF conducted emissions CISPR 11	Class B	The LCD is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

Electromagnetic immunity

The LCD Programming Base Station is intended for use in the electromagnetic environment specified below. The customer or the user of the system should ensure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment—guidance
Electrostatic discharge (ESD) IEC 61000-4-2	$\pm 2, 4, 6$ kV contact ± 8 kV air	N/A. The LCD does not have conductive surfaces. $\pm 2, 4, 6, 8$ kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	$\pm 2, 4, 6$ kV contact ± 8 kV air	± 2 kV for power supply lines ± 1 kV for input/output lines	Line power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	$\pm 0.5^*$, 1 kV differential mode ± 2 kV common mode	Line power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions, and voltage variations on power supply input lines IEC 61000-4-11	$< 5\% U_T$ ($> 95\%$ dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles $< 5\% U_T$ ($> 95\%$ dip in U_T for 5 sec)	$< 5\% U_T$ ($> 95\%$ dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles $< 5\% U_T$ ($> 95\%$ dip in U_T for 5 sec)	Line power quality should be that of a typical commercial or hospital environment. If the user requires continued operation during power line interruptions, it is recommended that the LCD be powered from the battery. Note U_T is the A.C. line voltage before application of the test level.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical commercial or hospital environment.
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms 10 KHz to 100 MHz*	Recommended separation distance: $d = 1.17 \sqrt{P}$

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment—guidance
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	10 V/m (compliance level adjusted to meet FDA limits) 26 MHz to 2.5 GHz Note At 80 MHz and 800 MHz, the higher frequency range applies.	Recommended separation distance: $d = 0.35 \sqrt{P}$ 80 MHz to 800MHz Recommended separation distance: $d = 0.70 \sqrt{P}$ 800MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 

Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the system is used exceeds the applicable RF compliance level above, the system should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the system.

Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Appendix: Limited warranty

Somnetics warrants the LCD to be free of defects in materials and workmanship and will perform in accordance with the product specifications for a period of 2 years from the date of sales by Somnetics to the dealer.

If the product fails to perform in accordance with the product specifications, Somnetics will repair or replace, at its option, any materials or parts of the LCD, which, upon Somnetics' examination appear defective. This does not cover damages caused by accident, misuse, abuse, alteration, and other defects not related to material or workmanship. Somnetics will pay customary freight charges from Somnetics to dealer location only.

Somnetics disclaims all liability for economic loss, loss of profits, overhead, or consequential damages which may be claimed to arise from any sale or use of this product. Some states do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation or exclusion may not apply to you.

This warranty is given in lieu of all other express warranties. In addition, any implied warranties, including warranty of merchantability or fitness for the particular purpose are limited to two years. Some states do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation or exclusion may not apply to you. This warranty gives you specific legal rights, and you may also have rights which vary from state to state.

To qualify for repair, replacement, or refund, the defective device must be returned to Somnetics within 30 days after the discovery of the defect. Any repair, replacement, or refund obligation would not apply if the device has been repaired or otherwise altered in a facility not authorized in writing by Somnetics To exercise your rights under this warranty, contact your local, authorized Somnetics dealer or Somnetics at 33 5th Avenue, New Brighton, Minnesota 55112 USA, 1.877.621.9626 or 1.651.621.1800.

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©2014 Somnetics International, Inc. All Rights Reserved.
ATLANTICO SYSTEMS LTD.