TRANSCEND[®] 365 mini **cpap**[™]

Quick Guide

For the
Transcend 365 miniCPAP Auto

Table of Contents

ENGLISH	3
FRANÇAIS	27
DEUTSCH	53
ITALIANO	80
ESPAÑOL	107
NEDERLANDS	133

Introduction

Thank you for choosing the Transcend 365 miniCPAP Auto™ by Transcend® - an obstructive sleep apnea medical device. The Transcend 365 miniCPAP Auto provides positive airway pressure in the range of 4-20 cmH₂O and includes pressure relief on exhalation with EZEX™ pressure relief technology.

- The Transcend 365 miniCPAP Auto's may be used every day on a bedside table and can be packed with you when you travel.
- An interactive LCD interface makes the 365 miniCPAP easy to operate.
- The Transcend 365 miniCPAP Auto is the only continuous positive airway pressure (CPAP) device licensed to use integrated Capillary Force Vaporization™ (CFV) humidification technology*. CFV instantaneously combines humidified vapor with the airstream when you inhale, providing a warm, moist, comfortable airflow. This humidification method requires less water than other humidification systems, which reduces the size of the humidifier.
- A DC power jack and USB port are also incorporated into the Transcend 365 miniCPAP Auto.

*Capillary Force Vaporization™ (CFV) humidification technology may be covered by one or more U.S. or foreign patents (patent 6,634,864, patent 7,431,570, patent 7,920,777, patent 8,201,752, patent 7,942,644, patent 5,008,207) and other U.S. or foreign patents pending.

Note: For more detailed descriptions, warnings, cautions and troubleshooting procedures refer to the Transcend 365 miniCPAP Auto User Manual that is available for download at www.MyTranscend.com/support or by calling Somnetics customer service at 877-621-9626 or info@somnetics.com to request a copy.

Components of the Transcend 365 miniCPAP Auto™

Begin by unpacking all items from the Transcend 365 miniCPAP Auto travel bag and inspecting them to ensure they were not damaged during shipment. Report any missing or damaged items to the home healthcare provider that sold the Transcend 365 miniCPAP Auto to you.

Included with Transcend 365 miniCPAP Auto

- Transcend 365 miniCPAP Auto device and integrated humidifier
- Universal power supply (PSA3) and power cords
- Travel bag
- Transcend 365 miniCPAP Auto Quick Guide
- USB cable

Accessories (sold separately)

- Air supply tube (compatible with standard 22 mm connector)
- P10 Battery™
- Patient Mask

Assembling the Transcend 365 miniCPAP Auto Device

Your mask connects directly to the Transcend 365 miniCPAP Auto via the air supply tube.

Follow these steps to assemble the Transcend 365 miniCPAP Auto:

1. Connect your mask (sold separately) to an air supply tube with a standard 22 mm connector (sold separately).
2. Connect the air supply tube directly to the Transcend 365 miniCPAP Auto device, making sure it is fully seated.
3. Plug the power supply into the Transcend 365 miniCPAP Auto. Place device on a flat, stable surface.
4. Connect the appropriate power cord to the power supply and connect the power cord to a wall outlet.

Setting up the Transcend 365 miniCPAP Auto

Home Screen

When supplied with power, you will see the LCD Home screen.


Center-Power, 'Select' button, Home, and Therapy Mode

Up- Settings

Down- Humidification, Ramp Activation


Left/Back- Compliance Report

Right- Last Session Report

The top right white numbers on the Home screen show the set pressure (large number) and the minimum and maximum pressures (small, lighter numbers - 4 and 20 in this image).

Control Panel


- Pressing one of the **four directional buttons** (Up, Down, Left or Right) on the Control Panel will immediately bring up the section you highlight on the Home screen. When in a section, these buttons help you navigate to highlight areas and make choices.
- Note that the **Left button is also the "back"** button. Press it to move back a screen or navigate back to the Home screen.
- Pressing the Center (hard plastic) button will start the PAP blower and send you to the Therapy Mode Screen. Press the down button to activate Ramp, if desired.

Note: When making choices in sections, navigate to highlight the value you want, then **press the Center button to save your choice.**

Using the LCD Screen


Therapy Settings

Enter the Settings screen by pushing the Up button to highlight the icon, then press the Center button to select it.

In the **Settings** section you will see four subsections:

- Therapy Settings
- Mask/Tube Settings
- Language/Units
- Device Information

A screenshot of the 'Patient Settings' menu. The menu is displayed in a list format with a dark header 'Patient Settings'. Below the header, four options are listed: 'Therapy Settings...', 'Mask/Tube Settings...', 'Language/Units...', and 'Device Information...'. The 'Therapy Settings...' option is highlighted with a grey background.


Press the Center button to enter the **Therapy Settings**.

Ramp Settings

The **Ramp** setting starts therapy pressure at a lower setting and gradually increases pressure until it reaches the prescribed therapy pressure, so you can fall asleep comfortably. To activate Ramp during therapy, press the down button.

From the Home screen, select the Settings icon using the directional buttons and the Center button. Toggle down to highlight Ramp and press center button to make comfort adjustments.

Use the Up or Down button to change the Ramp pressure (you cannot set it to a value less than the starting pressure). Remember to press the Center button to save your selection.


Use the right arrow to highlight and the Center button to select Ramp Duration in minutes. Press the Up or Down buttons to change the duration from 0 to 45 minutes, moving through increments of five minutes. **Press the Center button to save your selection.**

EZEX Settings

The EZEX feature allows you to reduce pressure when you exhale for more comfortable therapy.


Push the Down button to highlight EZEX, then press the Center button to adjust settings. The EZEX settings adjust the amount of pressure drop on exhalation, so that you can comfortably exhale while using your CPAP.

Use the Up or Down buttons to move the value from 0 to 3. A setting of 0 provides no exhalation relief, 3 provides the greatest exhalation relief. **Press the Center button to save your selection.**


Mask/Tube Settings

These settings help you select the appropriate mask and tube for optimal CPAP performance. Select the settings that apply for you. Use the Up/ Down buttons to highlight the setting, then press the Center button to make adjustments. Use the Up/ Down buttons to toggle through the setting options, then **press the Center button to save your selection.**

Language Settings

These settings help you select your native language and regional settings. Select the language that applies for you. **Press the Center button to save your selection.**

Use the Left/back button to toggle back through screen options or return to the Home screen.

Device Information


This screen shows device information.

Device Information			
SN	B18A20034	FW	2EC4
NE	1234	VN	2.0.0
BH	432.6	TH	405.8
HT	335.3	WH	687.8
HT◆	0.0	WH◆	0.0
FLT	0		
	000 000 000 000 000		
DT	2018/10/31 20:48:09		


Quick Reference Codes:

- SN- Device Serial Number
- FLT- Fault code summary. See Fault and Alert Code section in device user manual for more information.
- QR Barcode- link to current manuals and use resources.

Humidifier Settings

From the Home screen, push the Down button to highlight the Humidifier section. Push the Center button to choose it.

Use the Up or Down buttons to move the humidification level from 0 to 10 (0 = no humidification, 10 = maximum humidification). You will see the water droplet icons move to a higher or lower value. **Once you reach the desired setting, use the left and right buttons to toggle to the green check mark, if it's not already highlighted, then press the Center button to save your settings.** To exit out of the humidifier settings without saving any changes, select the red X by using the left and right buttons to highlight it, then press the Center button.


Push the Left/back button to return to the Home screen and you will see your new humidifier setting displayed on the Humidification icon.

Humidifier LCD indicators in Standby Mode

Water Status Indicators

These indicators will only be displayed if the humidifier is set to a level > 0. They will appear in lower right portion of the Home Screen

Flashing Blue Droplet - Out of Water


Flashing blue droplet

Indicated by a flashing blue droplet. Humidifier will be disabled during therapy. Can be caused by:

- No water in the reservoir – only displayed if heater/humidifier is set to level > 0.
- If this status appears but there is water in the reservoir, it indicates an error likely due to moisture in the CPAP. Allowing the system to dry out may correct the issue. To do this, run the CPAP in Dry Mode.

Note: After drying the system, disconnect and reconnect the power source to ensure the status is reset.

Flashing Yellow Droplet and Diamond – Tap Water Detected


Flashing yellow droplet

Only distilled water should be used with this humidifier. A flashing yellow droplet and yellow diamond indicate that water other than distilled water is detected in the reservoir. The humidifier can still be used, although doing so with any type of water other than distilled will shorten heater life and may void the warranty on the heater. This indicator may be caused by:

- Using water that is less pure than distilled water in the reservoir.
- Using distilled water in a reservoir that previously had tap water and wasn't sufficiently rinsed. If this is the case, it is recommended to dump out the reservoir water and rinse with distilled water. Replace water in reservoir with distilled water for use.
- Or, the water signal is corrupt and reporting an incorrect value. This is most likely caused by moisture in the CPAP. Allowing the system to dry out may correct the issue. To do this, run the CPAP in Dry Mode.

Note: After drying the system, disconnect and reconnect the power source to ensure the status is reset.

Flashing Red Droplet and X – Bad Water Signal Detected

Flashing red droplet and X

Indicated by a flashing red droplet and red 'X'. The humidifier will be disabled during therapy.


This indicates a corruption, most likely caused by moisture in the CPAP. Allowing the system to dry out may correct the issue. To do this, run the CPAP in Dry Mode.

Note: After drying the system, disconnect and reconnect the power source to ensure the status is reset.

Heater Status Indicators

These appear in lower right portion of the Home Screen (next to the humidification icon)

Flashing Red Exclamation (!) – Heater Over-Temperature Detected


Flashing red exclamation point

Indicated by a flashing red exclamation point.

- During the last therapy session, the heater reached a heat cutoff condition and turned off for the remainder of the session. The status will be reset when the next therapy session is started, and humidification will be enabled. If the status persists after subsequent therapy sessions, allowing the system to dry out may correct the issue. To do this, run the CPAP in Dry Mode.

Note: After drying the system, disconnect and reconnect the power source to ensure the status is reset.

Flashing Red X – No Heater Detected

Flashing red X

Indicated by a flashing red X. The humidifier will be disabled. Caused by:

- This indicates a corruption, most likely caused by moisture in the CPAP. Allowing the system to dry out may correct the issue. To do this, run the CPAP in Dry Mode.


Note: After drying the system, disconnect and reconnect the power source to ensure the status is reset.

If assistance is needed to reset an indicator, contact your medical equipment provider or Somnetics customer service 877-621-9626 or info@somnetics.com

Combination Indicators

The Water and Heater status symbols are not mutually exclusive. It is possible that both issues may occur and will be shown in combination in the lower right portion of the LCD screen.

Example:

Flashing red X


Humidifier LCD indicators in Therapy Mode

During therapy, the water/heater status and the humidifier setting is displayed in the lower right corner of the screen. The status consists of a water droplet that can have one of three colors:

- Gray droplet – Humidifier is disabled because humidifier setting is 0 (off).
- Blue droplet – Humidifier is enabled and operating.
- Red droplet – Humidifier is disabled because a condition was detected that does not permit humidifier operation.

Most users will be asleep and unaware if the status condition changes during the night. Additional status information is displayed on the Home screen in Standby (see section above) after therapy.


Therapy with Humidifier turned off (setting 0)

The gray droplet indicates that the humidifier is not active.


Therapy with Humidifier ON (setting > 0)

Normal operation, no errors, heater enabled.

The blue droplet indicates that the humidifier is active. Note: This includes cases where decreased water quality (tap water) was detected.

Error condition detected; humidifier disabled.

The red droplet indicates that the humidifier is not active. Normal CPAP therapy is not affected. This can be caused by any one of the following:

- Out of water condition detected.
- Bad water signal detected.
- Bad heater signal detected.
- Heater over-temperature detected.

More details about which error disabled the heater will be displayed on the standby screen (see LCD indicators in Standby Mode section) after therapy completes.

Standard user modes

Normal operation consists of four modes:

- Off** When the device is not connected to a power source, the device is off.
- Standby** When power is applied to the device the LCD screen lights up and the CPAP enters Standby Mode. Standby Mode is also initiated by pressing the power button when the device is in On Mode or if the mask is removed while in On Mode. As long as power is supplied to the device it will remain in Standby Mode until On or Drying Mode is initiated.
- On** When in On Mode, the blower is working and regulated device therapy pressure is being generated. The LED lights remain off. On Mode is initiated by pressing the center button when the device is in Standby Mode and the mask is worn by the patient.
- Drying Mode** This mode is accessed from Standby mode. Press and hold the Up and Down buttons at the same time for three seconds to start Drying Mode. This mode will last 30 minutes unless you interrupt it by pressing any button. Drying mode flushes your CPAP device, tubing and mask of residual moisture left from a therapy session.

Starting Therapy

Allow the Transcend 365 miniCPAP Auto to adjust to ambient temperature after removing it from storage. To begin therapy, plug in the power supply to your Transcend 365 miniCPAP Auto and press the center button to begin therapy. Press the down button to activate Ramp, if desired.

Note: To begin therapy from the home screen, first ensure that no icons are highlighted, then press the center button. You may also start therapy from any screen by holding the center power button for three seconds.

Press center button
to begin therapy.

Then press down
button to begin ramp.

Note: If the miniCPAP Auto™ device loses power while delivering therapy, the CPAP will power on as soon as power is restored. The user must hit the Center button to resume therapy.

Using the Ramp Function

1. Initiate Ramp by pressing the Down button after the Power button has been pressed (when CPAP is already delivering air pressure).
2. Pressure from the device will be lowered to a preset starting pressure and will gradually increase to your prescribed starting pressure.
3. To accelerate Ramp time, hold down the Down (Ramp) button. The device will quickly increase pressure until reaching your prescribed therapy pressure.

Ending therapy and entering Drying Mode

When you are finished with your therapy session, press the Center button to return your CPAP to Standby mode. Press and hold the Up and Down buttons for three seconds to start Drying Mode. This mode will last 30 minutes unless you interrupt it by pressing any button. After the 30-minute drying cycle, the blower will turn off and the device will automatically enter Standby Mode.

Dry mode flushes your CPAP device, tubing and mask of any residual moisture left from a therapy session.

Push and hold the Up and Down buttons to enter Drying Mode.

Caution: Do not attempt to use the Transcend 365 miniCPAP Auto while in Drying mode. The device will not provide the prescribed therapy pressure while in Drying mode.

The Transcend 365 miniCPAP Auto Humidifier

The Transcend 365 miniCPAP Auto is the only continuous positive airway pressure (CPAP) device licensed to use advanced, built-in Capillary Force Vaporization™ (CFV) humidification technology. CFV instantaneously combines humidified vapor with the airstream on inhalation for warm, moist, comfortable airflow. This unique humidifier utilizes a wick to draw distilled water up to the heating element, which then delivers warm vapor upon inhalation. This specialized method of delivering humidification requires less water, significantly reducing the size of the humidifier.

Caution! Use distilled water only. The use of distilled water is required for proper function of the heated humidifier. Failure to use distilled water will damage the humidifier and require repair or replacement. Alternatively, water may be filtered with the Transcend Portable Water Filter for use with your Transcend 365 miniCPAP humidifier. The ZeroWater™ Portable Filter Replacement Cartridge can be purchased directly from your Transcend dealer, or at most major retailers.

Filling the Reservoir

1. To remove the reservoir from the CPAP, first disconnect the power source to your CPAP, then set the device on a flat, stable surface and secure the reservoir with one hand while you lift and remove the CPAP with your other hand. **The pieces are**

simply joined by magnets- there is not a button or latch that releases the humidifier from the CPAP device.

2. Hold the reservoir at a slight angle and fill it with distilled water only. This step is best performed over a sink in case you spill water while filling the reservoir.

Caution: Do not perform this operation near electronic equipment. Always unplug the power source to your CPAP before attempting to remove the reservoir.


3. You can also use the Transcend 365 Water Filter Kit to fill your reservoir.

- Unpackage the Water Filter Kit and follow the instructions contained in the packaging to assemble the kit.
- Place the assembled kit on the Transcend 365 reservoir sitting on a flat surface.
- Add water to the fill line. Do not overflow the water filter as this will overflow the reservoir.
- Allow time for water to flow through the filter and fill the reservoir.

Note: It is essential to change the filter when indicated by the color change filter window (as shown below).

4. Check that the reservoir is at its full capacity by holding it vertically and viewing the water level. The water should be level with the max fill line.

Reservoir should be filled to the 'Max Fill' line.

Line up the heater wick with the opening in the reservoir.

Figure 1 Line the heater cartridge up with the opening in the reservoir. Lower the device onto the reservoir. Magnets will engage to hold the two parts together.

5. Once the reservoir is filled completely, set it on a secure, flat surface and line up the heater wick with the opening in the reservoir. Hold the reservoir in place and set the CPAP onto it so that the magnets snap the pieces together securely.

Note: Be careful while inserting the heater wick into the reservoir so that the reservoir seal does not get folded over, as this will generate water leaks.

Note: In the case that the reservoir seal becomes dislodged from the reservoir, call your Transcend miniCPAP dealer for a replacement.

Replacing the Heater Wick

Periodically, your heater wick will need replacement. If after using the device for several months you experience dry nose or throat symptoms, this may indicate the need to replace the wick. To order a replacement wick, call our customer service line at 877-621-9626.

1. Disconnect the power source, remove the humidifier reservoir, then turn the CPAP device over and you will see the wick cover. Turn it a half turn and pull to remove it.

2. Pull out the old wick to remove it and discard.
3. Place the new wick into the wick cover, then insert the wick into the humidifier.

Minimize handling of the new wick media to avoid contamination.

4. Lock the wick cover into place by giving it a half turn clockwise.

Heater Cartridge Replacement

If your heater cartridge should fail, you will see an alert notice on the LCD screen.

Note: To replace the heater cartridge, contact your CPAP dealer or Somnetics 877-621-9626 or info@somnetics.com for repair of your device.

Caring for your Transcend 365 miniCPAP Auto

The Transcend 365 miniCPAP Auto device is a maintenance-free device. However, it's important to keep it clean to ensure maximum efficiency and durability. Cleaning cycle frequencies identified below are checked against product service life.

Warning:

- Unplug the Transcend 365 miniCPAP Auto before cleaning.
- Do not submerge the Transcend 365 miniCPAP Auto or power supply in liquid.
- Prevent water from entering any openings of the device.
- Do not use harsh or abrasive cleaning agents to clean the device or any components.
- Do not attempt to sterilize the Transcend 365 miniCPAP Auto.
- Do not place cleaning materials, such as a cloth or liquid, into the device air inlet or air outlet connector.

Cleaning the CPAP Unit

Follow these instructions to clean the exterior of the Transcend 365 miniCPAP Auto.

Wipe the device with a dry, lint-free cotton cloth

Exterior cleaning of the device should be performed as follows:

Product	Periodic Cleaning Cycle	Product Service Life
Transcend 365 miniCPAP Auto	Weekly	5-Year

Caring for the Humidifier Reservoir


Follow these instructions to care for the humidifier reservoir.

1. Rinse the reservoir with distilled water weekly.
2. Reattach the humidifier reservoir to the PAP device.
3. Press and hold the Up and Down buttons for three seconds to start Drying Mode. This mode will last 30 minutes unless you interrupt it by pressing any button. After the 30-minute drying cycle, the blower will turn off and the device will automatically enter Standby Mode. Dry mode flushes your CPAP device, tubing and mask of any residual moisture.
4. Refill the reservoir with water prior to next use.

Replacing the Filter Media

During therapy, ambient air is drawn into the Transcend 365 miniCPAP Auto through an Air Inlet Filter. The Filter Media is very simple to replace. Simply turn the unit over and look for the door on one end that pops open easily.

Open door to remove the used filter media, replace with new filter media and secure the filter door. The filter media should be replaced after 6 months of use or sooner if it is dirty.

Questions? Contact Somnetics customer service 877-621-9626 or info@somnetics.com.

Appendix: Part Numbers

- Disposable parts
- Accessories
- Replacement parts

Note: All disposable parts, accessories and replacement parts are single patient reusable and are not for multi-patient use.

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Figure 2: Aligned the heating cartridge and the reservoir opening. Lower the device onto the reservoir. The magnets engage to hold the two parts together.

5. Once the reservoir is completely full, place it on a flat and safe surface and align the heating cartridge with the reservoir opening. Hold the reservoir in place and place the CPAP on top of it so that the magnets lock the pieces together securely.

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Somnetics® and Transcend are registered trademarks of Somnetics International, Inc. All rights reserved.
© 2019 Somnetics International, Inc.; all rights reserved.
Contact Somnetics International, Inc.
Somnetics International, Inc.
Atlantico Systems Ltd.