TRANSCEND[®] 365 mini **cpap**[™]

Quick Guide

Table of Contents

ENGLISH		3
FRANÇAIS		27
DEUTSCH		52
ITALIANO		77
ESPAÑOL		103
NEDERLANDS		128

Introduction

Thank you for choosing the Transcend 365 miniCPAP™ by Transcend® - an obstructive sleep apnea medical device. The Transcend 365 miniCPAP provides positive airway pressure in the range of 4-20 cmH₂O as prescribed by your clinician. The Transcend 365 miniCPAP product family includes the Transcend 365 miniCPAP and the Transcend 365 miniCPAP Auto.

- The Transcend 365 miniCPAP may be used every day on a bedside table and can be packed with you when you travel.
- An interactive LCD interface makes the Transcend 365 miniCPAP easy to operate.
- The Transcend 365 miniCPAP is the only continuous positive airway pressure (CPAP) device licensed to use integrated Capillary Force Vaporization™ (CFV) humidification technology*. CFV instantaneously combines humidified vapor with the airstream when you inhale, providing a warm, moist, comfortable airflow. This humidification method requires less water than other humidification systems, which reduces the size of the humidifier.
- A DC power jack and USB port are also incorporated into the Transcend 365 miniCPAP.

*Capillary Force Vaporization™ (CFV) humidification technology may be covered by one or more U.S. or foreign patents (patent 6,634,864, patent 7,431,570, patent 7,920,777, patent 8,201,752, patent 7,942,644, patent 5,008,207) and other U.S. or foreign patents pending.

NOTE: For more detailed descriptions, warnings, cautions and troubleshooting procedures refer to the Transcend 365 miniCPAP User Manual that is available for download at www.MyTranscend.com/support or by calling Somnetics customer service at 877-621-9626 or info@somnetics.com to request a copy.

Components of the Transcend 365 miniCPAP

Begin by unpacking all items from the Transcend 365 miniCPAP travel bag and inspecting them to ensure they were not damaged during shipment. Report any missing or damaged items to the home healthcare provider that sold the Transcend 365 miniCPAP to you.

Included with Transcend 365 miniCPAP

- Transcend 365 miniCPAP device and integrated humidifier
- Universal power supply (PSA3) and power cords
- Travel bag
- Transcend 365 miniCPAP Quick Guide
- USB cable

Accessories (Sold Separately)

- Air supply tube (compatible with standard 22 mm connector)
- P10 Battery™
- Patient Mask

Assembling the Transcend 365 miniCPAP Device

Your mask connects directly to the Transcend 365 miniCPAP via the air supply tube.

Follow these steps to assemble the Transcend 365 miniCPAP:

1. Connect your mask (sold separately) to an air supply tube with a standard 22 mm connector (sold separately).
2. Connect the air supply tube directly to the Transcend 365 miniCPAP device, making sure it is fully seated.
3. Plug the power supply into the Transcend 365 miniCPAP. Place device on a flat, stable surface.
4. Connect the appropriate power cord to the power supply and connect the power cord to a wall outlet.

Setting up the Transcend 365 miniCPAP

Home Screen

When supplied with power, you will see the LCD home screen.

Center - Power, 'Select' button, Home, and Therapy Mode

Up - Settings

Down - Humidification, Ramp


Left/Back - Compliance Report

Right - Last Session Report

The top right white numbers on the home screen show the set pressure (large number) and the minimum and maximum pressures (small, lighter numbers - 4 and 20 in this image).

Control Panel


- Pressing one of the **four directional buttons** (Up, Down, Left or Right) on the Control Panel will highlight the corresponding icon on the home screen. Pressing the Center button will then immediately bring up the highlighted section. When in a section, these buttons help you navigate to highlight areas and make choices.
- Note that the **Left button is also the "back"** button. Press it to move back a screen or navigate back to the home screen.
- Pressing the Center (hard plastic) button will start the PAP blower and send you to the Therapy Mode Screen.

NOTE: When making choices in sections, navigate to highlight the value you want, then press the Center button to save your choice.

Using the LCD Screen


Therapy Settings

Enter the Settings screen by pushing the Up button to highlight the icon, then press the Center button to select it.

In the **Settings** section you will see four subsections:

- Therapy Settings
- Mask/Tube Settings
- Language/Units
- Device Information


Press the Center button to enter the **Therapy Settings**.

Ramp Settings

The **Ramp** setting starts therapy pressure at a lower setting and gradually increases pressure until it reaches the prescribed therapy pressure, so you can fall asleep comfortably. Ramp will start automatically when therapy begins. To reactivate Ramp during the night, press the Down button.

From the home screen, select the Settings icon using the directional buttons and the Center button. Toggle down to highlight Ramp and press Center button to make comfort adjustments.

Use the Up or Down button to change the Ramp pressure (you cannot set it to a value less than the starting pressure). Remember to press the Center button to save your selection.


Use the Right button to highlight and the Center button to select Ramp Duration in minutes. Press the Up or Down buttons to change the duration from 0 to 45 minutes, moving through increments of five minutes. **Press the Center button to save your selection.**

EZEX Settings

The EZEX feature allows you to reduce pressure when you exhale for more comfortable therapy.


Push the Down button to highlight EZEX, then press the Center button to adjust settings. The EZEX settings adjust the amount of pressure drop on exhalation, so that you can comfortably exhale while using your CPAP.

Use the Up or Down buttons to move the value from 0 to 3. A setting of 0 provides no exhalation relief, 3 provides the greatest exhalation relief. **Press the Center button to save your selection.**


Mask/Tube Settings

These settings help you select the appropriate mask and tube for optimal CPAP performance. Select the settings that apply for you. Use the Up/Down buttons to highlight the setting, then press the Center button to make adjustments. Use the Up/Down buttons to toggle through the setting options, then **press the Center button to save your selection.**

Language Settings

These settings help you select your native language and regional settings. Select the language that applies for you. **Press the Center button to save your selection.** The Units option lets you change the pressure units from cmH₂O to hPa.

Use the Left/back button to toggle back through screen options or return to the home screen.

Device Information


This screen shows device information.

Device Information			
SN	B18A20034	FW	2EC4
NE	1234	VN	2.0.0
BH	432.6	TH	405.8
HT	335.3	WH	687.8
HT◆	0.0	WH◆	0.0
FLT	0		
	000 000 000 000 000		
DT	2018/10/31 20:48:09		


Quick Reference Codes:

- SN - Device Serial Number
- FLT - Fault code summary. See Fault and Alert Code section in device user manual for more information.
- QR Barcode - link to current manuals and use resources.

Humidifier Settings

From the home screen, push the Down button to highlight the Humidifier section. Push the Center button to choose it.

Use the Up or Down buttons to move the humidification level from 0 to 10 (0 = no humidification, 10 = maximum humidification). You will see the water droplet icons move to a higher or lower value. **Once you reach the desired setting, use the Left and Right buttons to toggle to the green check mark, if it's not already highlighted, then press the Center button to save your settings.** To exit out of the humidifier settings without saving any changes, select the red X by using the Left and Right buttons to highlight it, then press the Center button.

Note: If the device detects a detached or damaged humidifier it will display an error message in the Humidifier screen, indicating that the humidifier is not operational:

Humidifier LCD Indicators in Standby Mode

Water Status Indicators

These indicators will only be displayed if the humidifier is set to a level > 0. They will appear in lower right portion of the home Screen.

Flashing Blue Droplet - Out of Water

Flashing blue droplet

Indicated by a flashing blue droplet. Humidifier will be disabled during therapy. Can be caused by:

- No water in the reservoir – only displayed if heater/humidifier is set to level > 0.
- If this status appears but there is water in the reservoir, it indicates an error likely due to moisture in the CPAP. Allowing the system to dry out may correct the issue. To do this, run the CPAP in Dry Mode.

NOTE: After drying the system, disconnect and reconnect the power source to ensure the status is reset.

Flashing Yellow Droplet and Diamond – Tap Water Detected


Flashing yellow droplet and diamond

Only distilled water should be used with this humidifier. A flashing yellow droplet and yellow diamond indicate that water other than distilled water is detected in the reservoir. The humidifier can still be used, although doing so with any type of water other than distilled will shorten heater life and may void the warranty on the heater. This indicator may be caused by:

- Using water that is less pure than distilled water in the reservoir.
- Using distilled water in a reservoir that previously had tap water and wasn't sufficiently rinsed. If this is the case, it is recommended to dump out the reservoir water and rinse with distilled water. Replace water in reservoir with distilled water for use.
- Or, the water signal is corrupt and reporting an incorrect value. This is most likely caused by moisture in the CPAP. Allowing the system to dry out may correct the issue. To do this, run the CPAP in Dry Mode.

NOTE: After drying the system, disconnect and reconnect the power source to ensure the status is reset.

Flashing Red Droplet and X – Bad Water Signal Detected

Flashing red droplet and X


Indicated by a flashing red droplet and red 'X'. The humidifier will be disabled during therapy. This indicates a corruption, most likely caused by moisture in the CPAP. Allowing the system to dry out may correct the issue. To do this, run the CPAP in Dry Mode.

NOTE: After drying the system, disconnect and reconnect the power source to ensure the status is reset.

Heater Status Indicators

These appear in lower right portion of the home screen (next to the humidification icon)

Flashing Red X – Heater Detached or Damaged


Flashing red X

Contact your supplier for repair or replacement.

If assistance is needed to reset an indicator, contact your medical equipment provider or Somnetics customer service 877-621-9626 or info@somnetics.com

Combination Indicators

The Water and Heater status symbols are not mutually exclusive. It is possible that both issues may occur and will be shown in combination in the lower right portion of the LCD screen. Example:


Humidifier LCD Indicators in Therapy Mode

During therapy, the water/heater status and the humidifier setting is displayed in the lower right corner of the screen. The status consists of a water droplet that can have one of three colors:

- Gray droplet – Humidifier is disabled because humidifier setting is 0 (off).
- Blue droplet – Humidifier is enabled and operating.
- Red droplet – Humidifier is disabled because a condition was detected that does not permit humidifier operation.

Most users will be asleep and unaware if the status condition changes during the night. Additional status information is displayed on the home screen in Standby (see section above) after therapy.


Therapy with Humidifier Turned Off (Setting 0)

The gray droplet indicates that the humidifier is not active.

Therapy with Humidifier ON (Setting > 0)


Normal Operation, No Errors, Heater Enabled

The blue droplet indicates that the humidifier is active.

NOTE: This includes cases where decreased water quality (tap water) was detected.

Error Condition Detected; Humidifier Disabled

The red droplet indicates that the humidifier is not active. Normal CPAP therapy is not affected. This can be caused by any one of the following:

- Out of water condition detected.
- Bad water signal detected.
- Bad heater signal detected.
- Heater over-temperature detected.

More details about which error disabled the heater will be displayed on the standby screen (see LCD indicators in Standby Mode section) after therapy completes.

Standard User Modes

Normal operation consists of four modes:


Off	When the device is not connected to a power source, the device is off.
Standby	When power is applied to the device the LCD screen lights up and the CPAP enters Standby Mode. Standby Mode is also initiated by pressing the power button when the device is in On Mode or if the mask is removed while in On Mode. As long as power is supplied to the device it will remain in Standby Mode until On or Drying Mode is initiated.
On	When in On Mode, the blower is working and regulated device therapy pressure is being generated. On Mode is initiated by pressing the Center button when the device is in Standby Mode and the mask is worn by the patient.

Drying Mode This mode is accessed from Standby mode. Press and hold the Up button for three seconds to start Drying Mode. This mode will last 30 minutes unless you interrupt it by pressing the Center button to enter On mode or by pressing and holding the Up button for 3 seconds. Drying mode flushes your CPAP device, tubing and mask of residual moisture left from a therapy session.

Starting Therapy

Allow the Transcend 365 miniCPAP to adjust to ambient temperature after removing it from storage. To begin therapy, plug in the power supply to your Transcend 365 miniCPAP and press the Center button to begin therapy. If Ramp is configured, it will begin automatically. (Note: After therapy pressure has been reached, Ramp can be reactivated by pressing the Down button, if desired.)

NOTE: To begin therapy from the home screen, first ensure that no icons are highlighted, then press the Center button. You may also start therapy from any screen by holding the Center power button for three seconds.

NOTE: If the Transcend 365 miniCPAP device loses power while delivering therapy, the CPAP will power on as soon as power is restored and automatically resume therapy.

Transcend Desktop Software

The free Transcend Desktop Software can be downloaded at www.MyTranscend.com/support. The Transcend software allow you to view detailed compliance data and adjust the comfort features on your device.

TranSync

The TranSync Compliance Monitoring System is designed to provide patients with a secure and convenient way to view compliance data and transmit it to homecare providers and clinicians via the secure, cloud-based TranSync website (www.MyTranSync.com). Data can be uploaded to the TranSync website with the use of either the Transcend Desktop Software or the TranSync Bluetooth Module and mobile app.

Users must enter the email address associated with their TranSync account. If you don't already have a TranSync account, or have a profile set up by your homecare provider, visit www.MyTransync.com/Registration to set up your individual account before you attempt to transfer compliance information.

TranSync Mobile App

The TranSync mobile app (available for IOS and Android) can be downloaded from the Apple® and Google® app stores. The app user guides can be downloaded at www.MyTranscend.com/support. The app can be used to view compliance data uploaded to the TranSync cloud database (www.MyTranSync.com) via the desktop software or the TranSync Bluetooth Module.

TranSync Bluetooth Module

TranSync Mobile System Compatibility	
iOS	Version 11.1 or later
Android	Version 5.0 (API level 21) or later
Regulatory	FCC Part 15 Compliant 

To use your Bluetooth Module, insert with logo facing up.

To download the TranSync Mobile App or the TranSync Mobile App User Guide go to <https://MyTranscend.com/support>.

Using the Ramp Function

1. Pressure from the device will be lowered to a preset starting pressure and will gradually increase to your prescribed starting pressure.
2. To accelerate Ramp time, hold down the Down (Ramp) button. The device will quickly increase pressure until reaching your prescribed therapy pressure.

Ending Therapy and Entering Drying Mode

When you are finished with your therapy session, press the Center button to return your CPAP to Standby mode. Press and hold the Up button for three seconds to start Drying Mode. This mode will last 30 minutes unless you interrupt it by pressing the Center button to begin therapy or by pressing and holding the Up button for 3 seconds. After the 30-minute drying cycle, the blower will turn off and the device will automatically enter Standby Mode.

Dry mode flushes your CPAP device, tubing and mask of any residual moisture left from a therapy session.

Push and hold the Up button to enter Drying Mode

Caution: Do not attempt to use the Transcend 365 miniCPAP while in Drying mode. The device will not provide the prescribed therapy pressure while in Drying mode.

The Transcend 365 miniCPAP Humidifier

The Transcend 365 miniCPAP is the only continuous positive airway pressure (CPAP) device licensed to use advanced, built-in Capillary Force Vaporization™ (CFV) humidification technology. CFV instantaneously combines humidified vapor with the airstream on inhalation for warm, moist, comfortable airflow. This unique humidifier utilizes a wick to draw distilled water up to the heating element, which then delivers warm vapor upon inhalation. This specialized method of delivering humidification requires less water, significantly reducing the size of the humidifier.

Caution! Use distilled water only. The use of distilled water is required for proper function of the heated humidifier. Failure to use distilled water will damage the humidifier and require repair or replacement. Alternatively, water may be filtered with the Transcend Portable Water Filter for use with your Transcend 365 miniCPAP humidifier. The ZeroWater™ Portable Filter Replacement Cartridge can be purchased directly from your Transcend dealer, or at most major retailers.

Filling the Reservoir

1. To remove the reservoir from the CPAP, first disconnect the power source to your CPAP, then set the device on a flat, stable surface and secure the reservoir with one hand while you lift and remove the CPAP with your other hand. **The pieces are simply joined by magnets - there is not a button or latch that releases the humidifier from the CPAP device.**

2. Hold the reservoir at a slight angle and fill it with distilled water only. This step is best performed over a sink in case you spill water while filling the reservoir.

Caution: Do not perform this operation near electronic equipment. Always unplug the power source to your CPAP before attempting to remove the reservoir.


3. You can also use the Transcend 365 Water Filter Kit to fill your reservoir.

- Unpack the Water Filter Kit and follow the instructions contained in the packaging to assemble the kit.
- Place the assembled kit on the Transcend 365 reservoir sitting on a flat surface.
- Add water to the fill line. Do not overfill the water filter as this will overflow the reservoir.
- Allow time for water to flow through the filter and fill the reservoir.

NOTE: It is essential to change the filter when indicated by the color change filter window (as shown below).

4. Check that the reservoir is at its full capacity by holding it vertically and viewing the water level. The water should be level with the max fill line.

Reservoir should be filled to the 'Max Fill' line

5. Once the reservoir is filled completely, set it on a secure, flat surface and line up the heater wick with the opening in the reservoir. Hold the reservoir in place and set the CPAP onto it so that the magnets snap the pieces together securely.

NOTE: The humidifier reservoir must be emptied before transporting to avoid damage to the CPAP.


Figure 1 Line the heater cartridge up with the opening in the reservoir. Lower the device onto the reservoir. Magnets will engage to hold the two parts together.

NOTE: Be careful while inserting the heater wick into the reservoir so that the reservoir seal does not get folded over, as this will generate water leaks.

NOTE: In the case that the reservoir seal becomes dislodged from the reservoir, call your Transcend miniCPAP dealer for a replacement.

Replacing the Soft Heater Wick

The Transcend 365 miniCPAP will have either a ceramic wick or a replaceable soft wick included with the humidifier. The ceramic wick does not need to be replaced if distilled water is used. Should you have a replaceable soft wick and after using the device for several months you experience dry nose or throat symptoms, this may indicate the need to replace the wick. To order replacement wicks, call our customer service line at 877-621-9626.

1. Disconnect the power source, remove the humidifier reservoir, then turn the CPAP device over and you will see the wick cover. Turn it a half turn and pull to remove it.

2. Pull out the old wick to remove it and discard.
3. Place the new wick into the wick cover, then insert the wick into the humidifier.

NOTE: Minimize handling of the new wick media to avoid contamination.

4. Lock the wick cover into place by giving it a half turn clockwise.

Heater Cartridge Replacement

If your heater cartridge should fail, you will see an alert notice on the LCD screen.

NOTE: To replace the heater cartridge, contact your CPAP dealer or Somnetics 877-621-9626 or info@somnetics.com for repair of your device.

Caring for Your Transcend 365 miniCPAP

The Transcend 365 miniCPAP device is a maintenance-free device. However, it's important to keep it clean to ensure maximum efficiency and durability. Cleaning cycle frequencies identified below are checked against product service life.

Warning:

- Unplug the Transcend 365 miniCPAP before cleaning.
- Do not submerge the Transcend 365 miniCPAP or power supply in liquid.
- Prevent water from entering any openings of the device.
- Do not use harsh or abrasive cleaning agents to clean the device or any components.

- Do not attempt to sterilize the Transcend 365 miniCPAP.
- Do not place cleaning materials, such as a cloth or liquid, into the device air inlet or air outlet connector.

Cleaning the CPAP Unit

Follow these instructions to clean the exterior of the Transcend 365 miniCPAP.

Wipe the device with a dry, lint-free cotton cloth

Exterior cleaning of the device should be performed as follows:

Product	Periodic Cleaning Cycle	Product Service Life
Transcend 365 miniCPAP	Weekly	5-Year

Caring for the Humidifier Reservoir


Follow these instructions to care for the humidifier reservoir.

1. Rinse the reservoir with distilled water weekly.
2. Reattach the humidifier reservoir to the PAP device.
3. Press and hold the Up and Down buttons for three seconds to start Drying Mode.
This mode will last 30 minutes unless you interrupt it by pressing any button.
After the 30-minute drying cycle, the blower will turn off and the device will automatically enter Standby Mode. Dry mode flushes your CPAP device, tubing and mask of any residual moisture.
4. Refill the reservoir with water prior to next use.

Replacing the Filter Media

During therapy, ambient air is drawn into the Transcend 365 miniCPAP through an Air Inlet Filter. The Filter Media is very simple to replace. Simply turn the unit over and look for the door on one end that pops open easily.

Open door to remove the used filter media, replace with new filter media and secure the filter door. The filter media should be replaced after 6 months of use or sooner if it is dirty.

Questions? Contact Somnetics customer service 877-621-9626 or info@somnetics.com.

Appendix: Part Numbers

- Disposable parts
- Accessories
- Replacement parts

NOTE: All disposable parts, accessories and replacement parts are single patient reusable and are not for multi-patient use.

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Somnetics® and Transcend are registered trademarks of Somnetics International, Inc. All rights reserved.
© 2019 Somnetics International, Inc.; all rights reserved.
Contact Somnetics International, Inc.
Somnetics International, Inc.
Atlantico Systems Ltd.